




mono dosis


Skin Care Packaging



mono dosis
4 ml

Dropper Cap: HD-PE
Dropper: HD-PE
Bottle: PP


Plunger Cap: LD-PE



mono dosis
10 ml

Dropper Cap: HD-PE
Dropper: HD-PE
Bottle: PP

Plunger Cap: LD-PE



mono dosis
15 ml

Dropper Cap: HD-PE
Dropper: HD-PE
Bottle: PP

Plunger Cap: LD-PE

mono dosis



Press, Click and Open – die clevere Art der Applikation

Verpackungslösung für sensible Produkte – hygienisch und garantiert sicher verschlossen. Erst durch simples Drücken der Kappe öffnet sich das Packaging – das Produkt steht direkt zur Applikation bereit. Ein optimales Applikationssystem – wiederverschließbar und damit vielfältig einsetzbar. Ob für Einmalanwendungen, Ampullen-Produkte, Peeling, Puder oder Konzentrate.

Press, Click and Open – the clever way of application

Packaging solutions for delicate products – hygienic and a guaranteed safe closure. The packaging is initially opened with a simple pressing down of the cap – the product is ready for application. An optimal application system – resealable and therefore, repeatedly applicable. Whether for an ampoule-product, peeling, powder or a concentrate.